

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
ABO	ABO BLOOD TYPING	86900	N	-	Blood Bank	2 10ml pink	If baby draw pediatric tubes
ABORH	ABO AND RH	-	Y	-	Blood Bank	2 10ml pink	If baby draw pediatric tubes
<i>COMPONENT:</i>	BLOOD TYPING ABO	86900	-	1		2 10ml pink	If baby draw pediatric tubes
<i>COMPONENT:</i>	BLOOD TYPING RH D	86901	-	1		2 10ml pink	If baby draw pediatric tubes
ABSCRN	ANTIBODY SCREEN, RBC	86850	N	-	Blood Bank	2 10ml pink	If baby draw pediatric tubes
ACENTAB	ACENTROMERE B ABS	86235	N	-	LabCorp via MWMC	Serum	room temp
ACETA	ACETAMINOPHEN	82003	N	-	Chemistry	Serum	Drug Level, Do Not Use Gel Tube
ACHRAB	ACETYLCHO R BL AB	83520	N	-	LabCorp via MWMC	Serum	Refrigerate
ACKTQT	ACETON-KETON SERM QT	82010	N	-	Chemistry	Glass Red	
ACTH	ACTH ADRENOCORT	82024	N	-	Chemistry	Lavender	Spin specimen down
ACTNOMAB	ACTINOMYCES AB	86602	N	-	LabCorp via MWMC	serum	only stable 8 hours ambient, freeze or refrigerate
ADENOVIR	ADENOVIRUS AB	86603	N	-	LabCorp via MWMC	NP swab in media	BAL and CSF also ok
AFP	AFP SERUM	82105	N	-	LabCorp via MWMC	Serum	
AGCONV	ANGIOTENSIN CONV ENZ	82164	N	-	LabCorp via MWMC	Serum	refrigerate
AL	ALUMINUM	82108	N	-	LabCorp via MWMC	lavender	royal blue top and red tubes ok
ALA	ALA DELTA	82135	N	-			
ALB	ALBUMIN SERUM	82040	N	-	Chemistry	PST Green	
ALCOH	ALCOHOL (DRUG SCREEN SINGLE), URINE	G0431	N	-	Chemistry	Sterile Container	
ALDOLA	ALDOLASE	82085	N	-	LabCorp via MWMC	Serum	
ALDOS	ALDOSTERONE	82088	N	-	LabCorp via MWMC	Serum	
ALGAL	SPECTROPHOTO, NOS	84311	N	-			
ALLERGQT	ALLERGEN IGE QUANT	86003	N	-	LabCorp via MWMC	Serum	
ALP	ALKALINE PHOSPHATASE	84075	N	-	Chemistry	PST Green	
ALPISO	PHOSPHATASE, ALK ISO	84080	N	-	LabCorp via MWMC	serum	spin specimen down and separate as soon as blood clots
ALT	TRANSFERASE ALT	84460	N	-	Chemistry	PST Green	-
AMIOD	AMIODARONE QUANTITATIVE	80299	N	-	LabCorp via MWMC	serum	freeze
AMITRI	AMITRIPTYLINE	80152	N	-	LabCorp via MWMC	serum or plasma	Drug Level, cause for rejection GEL-barrier tubes
AMMBLD	AMMONIA	82140	N	-	Chemistry	PST Green	Put tube in ice
AMOXI	DRUG UNSPEC QT	80299	N	-			
AMUL	IMMUNOASSAY, OTHER	83516	N	-	LabCorp via MWMC		
AMY	AMYLASE	82150	N	-	Chemistry	PST Green	
ANA	ANTINUCLEAR AB (ANA)	86038	N	-	LabCorp via MWMC	serum	

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
ANACOMP	ANA COMPREHENSIVE PROFILE		Y	-	LabCorp via MWMC	serum	room temp
<i>COMPONENT:</i>	DNA AB, NATIVE-DOUBLE	86225	N	1			
<i>COMPONENT:</i>	ANTI-ENA	86235	N	8			
ANCAABTT	FLUORESCENT AB, TITER	86256	N	-			
ANCANEAB	ANCA NEUT CYTO AB		Y	-	LabCorp via MWMC	serum	refrigerate
<i>COMPONENT:</i>	FLUORESCENT AB, TITER	86256	-	3			
ANDRSTEN	ANDROSTENEDIONE	82157	N	-			
ANTENA	ANTI ENA	86235	N	-			
ANTNOMAB	ACTINOMYCES AB	86602	N	-	DUPLICATE		
ANTPLT	HEP INDC PLAT ANTI	86022	N	-	PHL	red or glass blue	sent stat 28060
ANTRNP	ANTI RNP	86235	N	-			
APHOAPH	ANTIPHOSPHATIDYLSERINE IGA/G/M (ANTIPHOSPHOLIPID A		Y	-	LabCorp via MWMC	serum	refrigerate
<i>COMPONENT:</i>	PHOSPHOLIPID AB	86148	-	3			
APTT	THROMBOPLAST TIMEPTT	85730	N	-	Coag	blue	
ARSENIC	ARSENIC, BLOOD	82175	N	-	LabCorp via MWMC	Royal blue top	send whole blood, room temp
ASOTITER	ASO, TITER	86060	N	-	LabCorp via MWMC	Serum	refrigerate
ASPERAG	ASPERGILLUS AG EIA (GALACTOMANNAN)	87305	N	-	PHL	Serum	
ASPERGAB	ASPERGILLUS AB	86606	N	-	LabCorp via MWMC	serum	refrigerate
ASPFUM	ASPERGILLUS FUMIGATU	86606	N	-	LabCorp via MWMC	serum	refrigerate
AST	TRANSFERASE AST	84450	N	-			
ASTXATX	ASSAY TOXIN/ANTITOXI	87230	N	-			
ATYPNEPN	ATYPICAL PNEUMO PNL	-	Y	-	PHL		
<i>COMPONENT:</i>	LEGIONELLA AB	86713	-	1			
<i>COMPONENT:</i>	MYCOPLASMA AB	86738	-	2			
AT3ACT	ANTITHROMBIN IIIACTI	85300	N	-			
AT3AG	ANTITHROMBIN III, AG	85301	N	-			
A1APHE	ALPHA 1-ANTITRYP PHE	82104	N	-	LabCorp via MWMC	Serum	refrigerate
A1APNL	ALPHA 1 ANTITRYP PNL	-	Y	-	LabCorp via MWMC		
<i>COMPONENT:</i>	ALPHA 1-ANTITRYP TOT	82103	-	1			
<i>COMPONENT:</i>	ALPHA 1-ANTITRYP PHE	82104	-	1			
A1AQT	ALPHA 1-ANTITRYP TOT	82103	N	-	LabCorp via MWMC	Serum	

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
BAGOTH	BACTERIAL AG-CSF-U-S	-	Y	-	LabCorp via MWMC	serum, csf, urine	2ml serum, 10 ml urine, 1ml csf
<i>COMPONENT:</i>	STREP B ASSAY W-OPTI	87802	-	1			
<i>COMPONENT:</i>	ALPHA 1-ANTITRYP PHE	87899	-	3			
BARTAB	BARTONELLA AB	86611	N	-	LabCorp via MWMC	serum	
BFCTWDI	CELL CT W-DIF BODY FLU	89051	N	-	Hematology	Lavender	
BFCTWODI	CELL CT W-O DIF BODY FL	89050	N	-	Hematology	Lavender	
BILD	BILIRUBIN DIRECT	82248	N	-	Chemistry	PST Green	
BILDNEO	BILIRUBIN NEO DIR	82248	N	-	Chemistry	PST Green	
BILITD	BILIRUBIN, TOT-DIR	-	Y	-			
<i>COMPONENT:</i>	BILIRUBIN TOTAL	82247	-	1			
<i>COMPONENT:</i>	BILIRUBIN DIRECT	82248	-	1			
BILT	BILIRUBIN TOTAL	82247	N	-	Chemistry	PST Green	
BILTNEO	BILIRUBIN NEO TOT	82247	N	-	Chemistry	PST Green	
BLDNOTFY	RBC BLOOD BANK 1 UNIT	-	-	-			
<i>PRIMARY CHARGE:</i>	S-H RED BLD-UNIT	P9021	N	-			
BMPWTCA	BASIC METABO CA TOTA	80048	N	-	Chemistry	PST Green	also red or sst acceptable
BNP	B-TYPE NATRIURETIC PEPTIDE	83880	N	-	Chemistry	Lavender	spin tube down
BORDPER	BORDETELL PERTUSS	87798	N	-	LabCorp via MWMC	Serum	refrigerate
BPRTAB	BORDETELLA AB	86615	N	-	LabCorp via MWMC	Serum	refrigerate
BRUCELAB	BRUCELLA AB	86622	N	-	LabCorp via MWMC	Serum	specimen should be free of bacterial contamination
BUN	UREA NITROGEN, QUAN (BUN)	84520	N	-			
B2GLYAB	BETA-2 GLYCOPROT AB	86146	N	-	LabCorp via MWMC	Serum	refrigerate
B2MIC	BETA-2 MICROGLOBULIN	82232	N	-	LabCorp via MWMC	Serum	refrigerate
B2TRS	IMM E-PHORS-URIN-CSF	86335	N	-			
CA	CALCIUM TOTAL	82310	N	-	Chemistry	PST Green	
CAD	CADMIUM	82300	N	-	LabCorp via MWMC	Royal blue top	whole blood, submit original tube
CAFR24HU	CATECHOL FRAC-24HRUR	82384	N	-	LabCorp via MWMC	30 mls urine	collect in brown urine collection container and add 30ml of 6HCl.
CAION	CALCIUM IONIZED	82330	N	-	Chemistry	Red	
CALCIT	CALCITONIN	82308	N	-	LabCorp via MWMC	Serum	Freeze
CALCULQT	CALCULUS ANAL QT	82360	N	-	LabCorp via MWMC	Calculi	Clean container or capsule. Specimen must be dry.
CALPRFEC	CALPROTECTIN, FECAL	83993	N	-	LabCorp via MWMC	Stool, Screw-cap plastic vial	refrigerate
CARBAMTO	CARBAMAZEPINE TOTAL	80156	N	-	Chemistry	Red or PST	avoid getl tubes. Tegretol

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
CARDGAM	CARDIOLIPIN AB G, A, M	-	Y	-			
<i>COMPONENT:</i>	CARDIOLIPIN AB, IGG	86147	-	3	PHL	Serum	27500 refrigerate
CARDIGM	CARDIOLIPIN IGM, AB	86147	N	-			
CARDLPG	CARDIOLIPIN AB, IGG	86147	N	-			
CATBLD	CATECHOLAMINES BLD	82383	N	-	LabCorp via MWMC	Plasma	lavender or heparin green Freeze
CA125	TUMOR AG CA 125	86304	N	-	LabCorp via MWMC	Serum	refrigerate
CA153	TUMOR AG CA 15-3	86300	N	-	LabCorp via MWMC	Serum	refrigerate
CA199	TUMOR AG CA 19-9	86301	N	-	LabCorp via MWMC	Serum	refrigerate
CA2729	CA 27-29	86300	N	-	LabCorp via MWMC	Serum	refrigerate
CBC	CBC W-PLT	85027	N	-	Hematology	Lavender	
CBCADIFF	CBC W-PLT AUTO COMPD	85025	N	-	Hematology	Lavender	
CBCMDIFF	CBC W-MANUAL DIFF	-	Y	-			
<i>COMPONENT:</i>	WBC DIFF, MAN	85007	-	1	Hematology	Lavender	
<i>COMPONENT:</i>	CBC W-PLT	85027	-	1	Hematology	Lavender	
CCP	CCP ANTIBODY	86200	N	-	LabCorp via MWMC	Serum	refrigerate or room temp
CDIFFAG	CLOSTRIDIUM DIF AG W/ REFLEX TOXIN IF POS	87449	N		PHL via MWMC		
CDIFFTX	CLOSTRIDIUM DIF TOXIN	87324	N	-	Microbiology	Stool	
CD4P	T CELL, ABSO CD4	86361	N	-	LabCorp via MWMC	whole blood	collect 7ml lavender and 7ml acda yellow top tube
CEA	CARCINOEMBRYONIC AG	82378	N	-	LabCorp via MWMC	Serum	refrigerate
CEDISAB	CELIAC DISEASE AB (PANEL)	-	Y	-			
<i>COMPONENT:</i>	MITOCHONDRIAL AB	86255	-	1			
<i>COMPONENT:</i>	GAMMAGLOB M-G-A-D	82784	-	1			
<i>COMPONENT:</i>	IMMUNOASSAY, OTHER	83516	-	4			
CERULO	CERULOPLASMIN	82390	N	-	LabCorp via MWMC	serum or plasma	refrigerate
CGA	CHROMOGRANIN A	86316	N	-	LabCorp via MWMC	serum or plasma	refrigerate
CHC	PCB/CHLORINATED PESTICIDES	82441	N	-			
CHEMASSA	CHEMILUM ASSAY NOS	82397	N	-			
CHLABT	CHLAMYDIA AB	86631	N	-	LabCorp via MWMC	Serum	refrigerate
CHLABTT	CHLAMYDIA AB TITTER	86631	N	-	LabCorp via MWMC	Serum	refrigerate
CHLAMYAB	CHLAMYDIA IGM AB	86632	N	-	LabCorp via MWMC	Serum	refrigerate
CHLGRNA	CHLAMYDIA-GC RNA	-	Y	-			
<i>COMPONENT:</i>	CHLMYD TRACH DNA AMP	87491	-	1			
<i>COMPONENT:</i>	N. GONORRHO, DNA, AMPPR	87591	-	1			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
CHLPSIGM	C. PSITTACI AB, IGM	86632	N	-			
CHLTDA	CHLMYD TRACH DNA AMP	87491	N	-			
CHLTDD	CHLMYD TRACH DNA DIR	87490	N	-			
CHOL	CHOLESTEROL TOTAL	82465	N	-	Chemistry	PST Green	
CHOLE	PSEUDOCHOLINESTERASE, SERUM	82480	N	-	LabCorp via MWMC	Whole blood	Lavender top tube refrigerate
CHORBC	CHOLINESTERASE RBC	82482	N	-	LabCorp via MWMC	Whole blood	Lavender top tube refrigerate
CHROM	CHROMIUM	82495	N	-	LabCorp via MWMC	Plasma	Plasma from royal blue top tube
CHROMQTM	CHROMATOG QUAN MULTI, URINE	82492	N	-			
CHROMTIS	CHROMOSOME ANALYSIS, TISSUE BIOPSIES (POC)	88262	N	-	LabCorp via MWMC	SEE WEBSITE	refrigerate
AH50	COMP, TOTAL CH50	86162	N	-	LabCorp via MWMC		
CK	CPK TOTAL	82550	N	-	Chemistry	PST, Red	
CKISO	CPK ISOENZYMES	82552	N	-			
CKISOPN	CPK ISOENZYMES	-	Y	-			
<i>COMPONENT:</i>	CPK TOTAL	82550	-	1			
<i>COMPONENT:</i>	CPK ISOENZYMES	82552	-	1			
CKMB	CPK MB FRACTION	82553	N	-	Chemistry	Red	
CL	CHLORIDE BLOOD	82435	N	-	Chemistry	PST Green	
CLINCHEN	CLIN CHEM TEST NOS (MISCELLANEOUS TEST)	84999	N	-			
CLSOTH	CHLORI OTHER FLD	82438	N	-			
ACT	COAG TIME ACTIVATED	85347	N	-			
CMP	COMP METAB PANEL	80053	N	-	Chemistry	PST Green	red and sst also acceptable
CMVDA	CYTOMEG DNA AMP PROB	87496	N	-			
CMVDQT	CYTOMEG DNA QUANT, PLASMA	87497	N	-	LabCorp via MWMC	Plasma from lav	Freeze
CMVIGG	CYTOMEGALOVIRUS AB, IGG (CMV	86644	N	-	LabCorp via MWMC	Serum	refrigerate
CMVIGM	CMV AB, IGM	86645	N	-	LabCorp via MWMC	Serum	refrigerate
COCAB	COCCIDIOIDES AB	86635	N	-	LabCorp via MWMC	Serum	refrigerate
COENZQ10	COENZYME Q10, TOTAL (CHROMATOG QUAN, S)	82491	N	-	LabCorp via MWMC	Serum	refrigerate
COLCHRQT	CHROMATOG QUAN SING (THIOPURINE METHYLTRANSF.)-T	82491	N	-			
COLDAGQT	COLD AGG, TITER	86157	N	-	LabCorp via MWMC	SST	keep tube in heat pack, put in incubator until clotted

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
COOXO2	CO-OXIMETRY-O2 SAT	-	Y	-	RT-ABG		
<i>COMPONENT:</i>	CARBON MONOXIDE QUAN	82375	-	1			
<i>COMPONENT:</i>	BLD GAS W-O2 SAT	82805	-	1			
<i>COMPONENT:</i>	METHEMOGLOBIN, QUAN	83050	-	1			
<i>COMPONENT:</i>	HEMOGLOBIN	85018	-	1			
COPPER	COPPER SERUM	82525	N	-	LabCorp via MWMC	Serum or plasma	serum from red or plasma from royal blue top tube
CORTAM	CORTISOL, SERUM AM	82533	N	-	Chemistry	red	
CORTPM	CORTISOL SERUM, PM	82533	N	-	Chemistry	red	
CORTT	CORTISOL TOTAL	82533	N	-	Chemistry	red	
COXSAA	COXSACKIE A AB	86658	N	-	LabCorp via MWMC	serum or plasma	red, lavender, heparin pst or blue tops all ok
CO2	CARBOX DIOXIDE (CO2)	82374	N	-			
CPEP	C-PEPTIDE	84681	N	-	LabCorp via MWMC	serum	refrigerate
CRAYG	CRYPTOCOCCAL AG SER (AGENT NOS ASSAY)	87899	N	-	LabCorp via MWMC	2ml serum	refrigerate
CREA	CREATININE BLOOD	82565	N	-	Chemistry	PST green	
CREATCLR	CREATININE CLEARANCE	82575	N	-	Chemistry	blood and urine	
CREATINE	CREATINE	82540	N	-	LabCorp via MWMC	Serum	Freeze
CREOTH	CREATININE, FLUID	82570	N	-	Chemistry	PST Red	sterile container ok too
CRP	CRP	86140	N	-	Chemistry	Red	
CRPHS	CRP, HIGHLY SENS	86141	N	-	LabCorp via MWMC	Serum or plasma	refrigerate
CRYAG	ORGANISM IDENTIFICATION, YEAST	87106	N	-			
CRYAGEIA	CRYPTOSPORIDIUM AG EIA	87328	N	-	Microbiology	Stool	
CRYOF	CRYOFIBRINOGEN	82585	N	-	LabCorp via MWMC	Plasma	blue, lavender, or grey top tubes
CRYOG	CRYOGLOBULIN	82595	N	-	LabCorp via MWMC	3 ml Serum	refrigerate
CRYONTFY	CRYOPRECIPITATE 1 UNIT	-	-	-	Blood Bank		
<i>PRIMARY CHARGE</i>	S-H CRYO EACH UNIT	P9012	N	-			
CRYPSMR	CRYPTOSPORIDIUM SMEAR	87206	N	-			
CRYPTEIA	CRYPTOSPORIDIUM EIA, STOOL	87328	N	-	Microbiology		
CRYPTOAG	FUNGI IDENTIF YEAST	87106	N	-			
CSFALB	ALBUMIN CSF	82042	N	-	Chemistry	green or red	
CSFCRYAG	CRYPTOCOCCAL AG CSF (PART AGGLUT.SCREEN)	86403	N	-			
CSFCTDIF	CSF CELL CT W-DIFF	89051	N	-	Hematology	lavender	
CSFGLU	GLUCOSE CSF	82945	N	-	Chemistry	green or red	
CSFMYEPR	MYELIN BASIC PROTEIN	83873	N	-	LabCorp via MWMC	1ml CSF	FREEZE

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
CSFPROQT	PROTEIN CSF QT	84157	N	-	Chemistry	sterile container	
CSFPRTTEL	PROTEIN TOT ELECTROPHORESIS CSF	-	Y	-			
<i>COMPONENT:</i>	PROT E-PHOR-URIN-CSF	84166	-	1			
<i>COMPONENT:</i>	PROTEIN TOT OTHR SRC	84157	-	1			
CSFVDRL	VDRL, CSF	86592	N	-	PHL	.5ml CSF	FREEZE send stat
CULAE	CULTURE BACTERIA OTH	-	-	-			
<i>PRIMARY CHARGE:</i>	CULTURE BACTERIA OTH	87070	N	-			
CULAFBB	CULTURE BLOOD, AFB (#64500)	87116	N	-	PHL		CALL FOR BLOOD CULTURE BOTTLES
CULAFBPN	CULTURE AFB-DIR SMEAR	-	Y	-			
<i>COMPONENT:</i>	CULT MYCOBACTERIA	87116	-	1			
<i>COMPONENT:</i>	SMEAR FLUORES-ACI ST	87206	-	1			
CULAN	CULTURE BACTERIA ANA ANY SO	87075	N	-			
<i>PRIMARY CHARGE:</i>	CULTURE BACTERIA ANA ANY SO	-	-	-			
CULBLFNG	CULTURE FUNGUS BLOOD	87103	N	-			
CULCATP	CULTURE CATHETER TIP	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT CATHETER TIP	87070	N	-			
CULCMVTC	CULTURE CMV TISSUE CULT (VIRUS INOC.)	87252	N	-			
CULCSF	CULTURE CSF	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT CSF	87070	N	-			
CULEAR	CULTURE EAR	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT EAR	87070	N	-			
CULECOL	CULTURE E COLI	-	-	-			
<i>PRIMARY CHARGE:</i>	CULTURE E COLI	87046	N	-			
CULFNGBF	FUNGI IDENTIF, MOLD	87107	N	-			
CULFNGIS	CULTURE FUNGUS ISOLATION	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT FUNGUS ISOLATIO	87102	N	-			
CULGC	CULTURE GC	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT GC	87070	N	-			
CULMRSA	CULTURE SCREEN MRSA	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT SCREEN MRSA	87081	N	-			
CULRSP	CULTURE RESPIRATORY	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT RESPIRATORY	87070	N	-			
CULSTB	CULTURE GROUP B STREP	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT GROUP B STREP	87081	N	-			
CULSTL	CULTURE STOOL BACT W-ISOL	-	-	-			
<i>PRIMARY CHARGE:</i>	CULT STL BACT W-ISOL	87045	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
CULSTLAD	CULT STL BACT EA ADD	87046	N	-			
CULTBA	CULTURE BRONCH ASPIRATE	-	-	-			
PRIMARY CHARGE:	CULTURE BRONCH ASPIRATE	87070	N	-			
CULTBF	CULTURE BODY FLUID	-	-	-			
PRIMARY CHARGE:	CULTURE BODY FLUID	87070	N	-			
CULTBLD	CULTURE BLOOD	-	-	-			
PRIMARY CHARGE:	CULTURE BLOOD BACTER	87040	N	-			
CULTCAMP	CULTURE CAMPYLOBACTER	-	-	-			
PRIMARY CHARGE:	CULT CAMPYLOBACTER	87046	N	-			
CULTCHL	CULTURE CHLAMYDIA	87110	N	-			
CULTENTV	VIRUS ISOLATION (VIRUS INOCULATION TISS)	87252	N	-			
CULTEYE	CULTURE EYE	-	-	-			
PRIMARY CHARGE:	CULTURE EYE	87070	N	-			
CULTGAS	CULTURE GRP A STREP	-	-	-			
PRIMARY CHARGE:	CULT GRP A STREP	87081	N	-			
CULTGEN	CULTURE GENITAL	-	-	-			
PRIMARY CHARGE:	CULT GENITAL	87070	N	-			
CULTIMG	CULTURE TYPE, IMMUNOLOGI	87147	N	-			
CULTIS	TISSUE HOMOGENIZ CUL	87176	N	-			
CULTLEGI	LEGIONELLA CULTURE	87081	N	-			
CULTNAS	CULTURE NASAL	-	-	-			
PRIMARY CHARGE:	CULT NASAL	87070	N	-			
CULTNUCA	CULT TYPE, NUCLEIC AC (AFB ID)	87149	N	-			
CULTTHRO	CULTURE THROAT	-	-	-			
PRIMARY CHARGE:	CULT THROAT	87070	N	-			
CULTTIBI	CULTURE TISSUE BIOPSY	-	-	-			
PRIMARY CHARGE:	CULT TISSUE BIOPSY	87070	N	-			
CULTUR	CULTURE URINE COL CT	-	-	-			
PRIMARY CHARGE:	CULT URINE-COLONY CT	87086	N	-			
CULTWND	CULTURE WOUND	-	-	-			
PRIMARY CHARGE:	CULTURE WOUND	87070	N	-			
CULYERS	CULTURE YERSENIA	-	-	-			
PRIMARY CHARGE:	CULTURE YERSENIA	87046	N	-			
CYANIDE	CYANIDE	82600	N	-	PHL	Whole blood	10 heparin green or 4ml lav. Do not spin. 82940
CYCLOS	CYCLOSPORINE	80158	N	-	LabCorp via MWMC	Whole blood	lavender or sodium heparin green
CYTOPCR	CYTOMEG PCR, BLOOD	87496	N	-	LabCorp via MWMC	Whole blood	4ml Lav, or yellow acd tube spun down and send plasma

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
C1	C1 COMPLEMENT	86160	N	-	LabCorp via MWMC	Serum	separate from cells within an hour and keep at room temp
C1ESIF	COMPLEMENT C1 ESTERASE INHIBITOR, FUNCTIONAL	86161	N	-	LabCorp via MWMC	Serum	Freeze
C1EST	C1 ESTERSE INHIBITOR, SERUM	86160	N	-	LabCorp via MWMC	Serum	separate from cells within an hour and keep at room temp
C3	C3 COMPLEMENT	86160	N	-	LabCorp via MWMC	Serum	remove from cells with in 45 mins and refrigerate
C4	C4 COMPLEMENT	86160	N	-	LabCorp via MWMC	Serum	allow to clot for 15-30 mins before separating, refrigerate
DAMULDNA	BORDETELLA PERTUSS DET AGT D	87801	N	-	PHL	nasopharyngeal	58466 place swab into universal transport media and refrigerate
DANOSDNQ	DETECT AGENT, DNAQ	87799	N	-			
DANSDNA	DETECT AGENT NOS/DNA	87797	N	-			
DANSDNAM	NOROVIRUS PCR (DETECT AGT NOS)	87798	N	-	LabCorp via MWMC	stool	refrigerate
DAT	DIRECT COOMBS TEST	86880	N	-	Blood Bank	whole blood	
DDIMQT	D-DIMER, QT	85379	N	-	Coag	Blue top	
DHEA	DHEA	82626	N	-	LabCorp via MWMC	Serum	refrigerate
DHEAS	DHEA-S	82627	N	-	LabCorp via MWMC	Serum	refrigerate
DIGITOX	DIGITOXIN (DRUG UNSPEC QUANT)	80299	N	-			
DIGOXIN	DIGOXIN	80162	N	-	Chemistry	Serum	
DIHYDTES	DIHYDROTTESTOSTERONE	82651	N	-	LabCorp via MWMC	Serum	Freeze
DIHYDX25	VIT D 1, 25 DIHYDROX	82652	N	-	LabCorp via MWMC	Serum	Protect from light
DILFRTOT	PHENYTOIN (DILANTIN) FREE-TOTAL	-	Y	-			
COMPONENT:	PHENYTOIN TOTAL	80185	-	1			
COMPONENT:	PHENYTOIN FREE	80186	-	1			
DNAAB	DNA AB, NATIVE/DOUBL	86225	N	-			
DRSCINFA	DRUG SCRIN, QUAL-MULTI	-	N	-			
DRSCR10	DRUG SCREEN 10	-	Y	-	Chemistry	urine	
COMPONENT:	UDS CHEM ANALYZER	G0431	-	10			
DRSCR5	DRUG SCREEN 5	-	-	-			
COMPONENT:	UDS CHEM ANALYZER	G0431	Y	5			
DRSQLMC	UDS MULTI QL TRIAGE	G0431	N	-			
EBCSTLOU	ENCEPHALITIS, ST LOUIS	86653	N	-	LabCorp via MWMC	Serum	room temp
EBPN1	EPSTEIN-BARR AB, EA	86663	N	-	LabCorp via MWMC	Serum	room temp
EBVIGGVC	EBV IGG, VCA	86665	N	-	LabCorp via MWMC	Serum	room temp

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
EBVIGMVC	EBV IGM, VCA	86665	N	-	LabCorp via MWMC	Serum	room temp
EBVPCRQL	EPSTEIN-BARR VIRUS PCR, QL	86665	N	-	LabCorp via MWMC	whole blood	whole blood from 4ml lav or csf room temp
ECHINOAB	ECHINOCOCCUS AB	86682	N	-	LabCorp via MWMC	Serum	refrigerate or freeze
EHISTAB	E. HISTOLYTICA AB (AMEBIASIS)	86753	N	-	LabCorp via MWMC	Serum	refrigerate
ENCCA	ENCEPHALITIS, CALIF	86651	N	-	LabCorp via MWMC	Serum	room temp
ENCEEQ	ENCEPHALITIS, EAST EQUINE	86652	N	-	LabCorp via MWMC	Serum	room temp
ENCEWEQ	ENCEPHALITIS, WEST EQUINE	86654	N	-	LabCorp via MWMC	Serum	room temp
ENCSTLOU	ENCEPHALITIS, ST LOUIS	86653	N	-	LabCorp via MWMC	Serum	room temp
ENDOMAB	ENDOMYSIAL IGA AB	86255	N	-	LabCorp via MWMC	Serum	room temp
ENTERPCR	ENTEROVIRUS RT-PCR	87498	N	-	LabCorp via MWMC	CSF, NP, SWAB	refrigerate
ENTVAB	ENTEROVIRUS AB	86658	N	-			
EOSNASSM	EOS SMEAR	89190	N	-			
ERYTHR	ERYTHROPOIETIN	82668	N	-	LabCorp via MWMC	Serum	refrigerate
ESOBDFL	EOSINOPHILS BODY FLUID	89050	N	-			
ESRM	SED RATE NONAUTOMATE	85651	N	-	Hematology	Lavender	
ESTRA	ESTRADIOL	82670	N	-	LabCorp via MWMC	Serum	refrigerate
ESTROGFR	ESTROGENS FRACT	82671	N	-	PHL	Serum	41621 refrigerate
ETHANOL	ALCOHOL, SERUM	82055	N	-	Chemistry	red or green	
ETHGLYC	ETHYLENE GLYCOL	82693	N	-	LabCorp via MWMC	Serum	original tube either sst or heparin green, refrigerate
EUGLBY	EUGLOBULIN LYSIS	85360	N	-	PHL 27100	4.5 glass blue	phl must receive within 30 mins. Or send 2ml plasma frozen
FACIX	FACTOR IX	85250	N	-	PHL	4.5 glass blue	27320 send stat
FACTV	FACTOR V ASSAY	85220	N	-	PHL	4.5 glass blue	27280 send stat
FACVIII	FAC VIII, VON WILLEBRAND				PHL	4.5 glass blue	27240 send stat
FACVIIN	FACT VIII, VON WILLEBRAND, MULTIMERS	85247	N	-	PHL	4.5 glass blue	97165 send 2ml plasma frozen stat
FACXI	FACTOR XI	85270	N	-	PHL	4.5 glass blue	27360 send stat room temp
FACXII	FACTOR XII	85280	N	-	PHL	4.5 glass blue	27400 send stat room temp
FAC5LPC	FACTOR V LEIDEN PCR	-	Y	-	PHL via MWMC	4.5mL lavender	refrigerate
COMPONENT:	MOLEC ISO EXT HP EA	83891	-	1			
COMPONENT:	MOLEC NUC ACID PROBE	83896	-	2			
COMPONENT:	MOLEC NUC AMPLI EA	83898	-	1			
COMPONENT:	MOLEC MUTATION SCAN	83903	-	1			
COMPONENT:	MOLEC INTERP-RPT	83912	-	1			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
FAC8RE	FACTOR VIII REL AG	85244	N	-			
FAC8RIST	FACTOR VIII RISTOCET	85245	N	-			
FCFTLPQL	FAT-LIPID FECES QUAL	82705	N	-			
FCFTLPQT	FAT-LIPID FECES QUAN	82710	N	-			
FCGIA	GIARDIA AG STOOL	87329	N	-	Microbiology		
FCLACTQL	LACTOFERRIN, FEC QAL	83630	N	-			
FCLYTES	ELECTROLYTES, FECES	-	Y	-			
<i>COMPONENT:</i>	CHLORIDE OTHER FLUID	82438	-	1			
<i>COMPONENT:</i>	CLIN CHEM TEST NOS	84999	-	1			
<i>COMPONENT:</i>	SODIUM OTHER SOURCE	84302	-	1			
FDPAGGSL	FDP-FSP AGG SLIDE	85362	N	-			
FECFAT	FAT STAIN FECAL RD	89125	N	-			
FENA	FRACTIONAL EXCRETION OF SODIUM	-	Y	-	Chemistry		
<i>COMPONENT:</i>	CREATININE BLOOD	82565	-	1			
<i>COMPONENT:</i>	CREATININE UR-OTHER	82570	-	1			
<i>COMPONENT:</i>	SODIUM, SER PL OR WB	84295	-	1			
<i>COMPONENT:</i>	SODIUM, URINE RANDOM	84300	-	1			
FERRIT	FERRITIN	82728	N	-	Chemistry		
FETFIBRO	FETAL FIBRONECTIN	82731	N	-			
FETHGBK	FETAL HGB KLEIHAUER-BETKE	85460	N	-			
FETHGBRS	FETAL HGB, ROSETTE	85461	N	-			
FETIBC	TOT IRON BINDING	-	Y	-	Chemistry		
<i>COMPONENT:</i>	IRON	83540	-	1			
<i>COMPONENT:</i>	IRON BINDING CAPACTY	83550	-	1			
FETLUNGR	FETAL LUNG, L-S RATIO	83661	N	-			
FFPNOTFY	FFP BLOOD BANK 1 UNIT	P9017	N	-			
FIB	FIBRINOGEN	85384	N	-			
FLDALB	ALBUMIN FLUID	82042	N	-	Chemistry	red or green	sterile container also ok. Urine sent out as UALB
FLDAMY	AMYLASE FLUID	82150	N	-	Chemistry	red or green	sterile container also ok. Avoid pst and sst
FLDAST	TRANSFERASE AST, BODY FLUID	84450	N	-			
FLDBIL	BILIRUBIN FLUID	82247	N	-			
FLDCEA	CARCINOEMBRYONIC AG, FLUID	82378					
FLDCHOL	CHOLESTEROL FLUID	82465	N	-			
FLDGLU	GLUCOSE FLUID	82945	N	-			
FLDLA	LACTIC ACID, BODY FLUID	83605	N	-			
FLDLDH	LDH FLUID	83615	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
FLDLIP	LIPASE FLUID	83690	N	-			
FLDOSML	OSMOLALITY FLUID	83930	N	-			
FLDPH	PH FLUID	83986	N	-			
FLDRFQT	RHEUMATOID FACT QUAN	86431	N	-			
FLDTRIG	TRIGLYCERI FLUID	84478	N	-			
FLDURIC	URIC ACID FLUID	84560	N	-			
FLECAI	FLECAINIDE TAMECORR	80299	N	-			
FLFPIA	FETAL LUNG MATURITY CASCADE	83663	N	-			
FLOABSCR	FLUORESCENT AB, SCREE	86255	N	-			
FLOABTIT	FLUORESCENT AB, TITER	86256	N	-			
FLUABAB	INFLUENZA A-B AB SER	-	Y	-			
<i>COMPONENT:</i>	INFLUENZA VIRUS AB	86710	-	2			
FLUABAG	INFLUENZA A-B, AG, EIA	87400	N	-			
FLUH1N1	FLU A NOVEL H1N1-PCR	87798	N	-			
FOBT	FECED OCCULT BLOOD 1-3	82270	N	-			
FOBTOTH	FECED OCCULT BLOOD, GASTRIC	82271	N	-			
FOLRCHCT	FOLATE, RBC W-HCT	-	Y	-			
<i>COMPONENT:</i>	FOLIC ACID RBC	82747	-	1			
<i>COMPONENT:</i>	HEMATOCRIT, OTHER	85014	-	1			
FOLSER	FOLIC ACID SERUM	82746	N	-			
FTULARAB	FRANCISELLA TULARENS-IGG & IGM	86668	N	-			
FRUCT	GLYCATED PROTEIN (FRUCTOSAMINE)	82985	N	-			
FSH	FOLLICLE STIM HORM	83001	N	-			
FTAABIF	FTA ABS IFA SERMIG	86780	N	-			
FT3	T3, FREE	84481	N	-			
FT4	THYROXINE T4, FREE	84439	N	-			
AGDNOSEM	DETEC NOS, EIA MUL (FUNGITELL)	87449	N	-	LabCorp via MWMC	Serum	http://www.viracor.com/tmp/Fungitell-Assay-Sheet.pdf
GAD65AB	GAD-65 AB	83519	N	-			
GASTRI	GASTRIN	82941	N	-			
GBMAB	GBM AB	83516	N	-			
GBMABIG	BASMENT MB IGG GBMAB	83520	N	-			
GCDNAQL	N. GONORRHO, DNA, AMPPR	87591	N	-			
GENTPK	GENTAMICIN PEAK	80170	N	-			
GENTR	GENTAMICIN RANDOM	80170	N	-			
GENTTR	GENTAMICIN TROUGH	80170	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
GGT	GGT	82977	N	-			
GHAB	GROWTH HORMONE AB	86277	N	-			
GLADIGG	GLIADIN AB IGG	83516	N	-			
GLIADIGA	GLIADIN AB IGA	83516	N	-			
GLU	GLUCOSE, QT	82947	N	-			
GTT-2H	GLUCOSE TOLERANCE 2HR - ADA	-	Y	-			
<i>COMPONENT:</i>	GLUCOSE TOL-FASTING FOR GTT2HR	82947	-	1			
<i>COMPONENT:</i>	GLUCOSE TOL-2HR FOR GTT2HR	82947	-	1			
GLUTOL2H	GLUCOSE TOLERANCE 2HR - GESTATIONAL SCREEN	-	Y	-			
<i>COMPONENT:</i>	GLUCOSE TOL-FASTING FOR GTT2HR		-	-			
<i>COMPONENT:</i>	GLUCOSE TOL-1HR FOR GTT2HR	82951	-	1			
<i>COMPONENT:</i>	GLUCOSE TOL-2HR FOR GTT2HR		-	-			
GLUTOL3H	GLUCOSE TOLERANCE 3HR - GESTATIONAL SCREEN	-	Y	-			
<i>COMPONENT:</i>	GLUCOSE TOL-FASTING FOR GTT3HR	82952	-	1			
<i>COMPONENT:</i>	GLUCOSE TOL-1HR FOR GTT3HR	82951	-	1			
<i>COMPONENT:</i>	GLUCOSE TOL-2HR FOR GTT3HR	-	-	-			
<i>COMPONENT:</i>	GLUCOSE TOL-3HR FOR GTT3HR	-	-	-			
GRAMST	SMEAR, GRAM STAIN	87205	N	-			
GRH	GROWTH HORMONE, SERUM (EA SPEC)	87205	N	-			
GTDNAHC	GENOTYPE, DNA, HEPATITIS C - NONREFLEX	87902	N	-			
GTT1HPD	GLUCOSE, POST DOSE	82950	N	-			
G6PD	G6PD	-	Y	-			
<i>COMPONENT:</i>	G6PD, QUANT	82955	-	1			
<i>COMPONENT:</i>	RBC COUNT	85041	-	1			
G6PDQT	G6PD, QT	82955	N	-			
HACUTEPN	HEPATITIS ACUTE PANEL	80074	N	-			
HAIGM	HEPATITIS A AB IGM	86709	N	-			
HAPF	HEPATITIS A AB, TOTAL	86708	N	-			
HAPTOQT	HAPTOGLOBIN, QT	83010	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
HBCABT	HEP B CORE AB, TOTAL	86704	N	-			
HBCIGM	HEPATITIS B CORE AB IGM	86705	N	-			
HBDNAA	HEPATITIS B, DNA AMPP	87516	N	-			
HBDNAQT	HEPATITIS B DNA DETECT QT	87517	N	-			
HBEAB	HEPATITIS BE AB (HBEAB)	86707	N	-			
HBEAG	HEPATITIS B ENVELOPE AG	87350	N	-			
HBSAB	HEP B SURFACE AB	86706	N	-			
HBSAG	HEPATITIS B SURFACE AG (HBSAG)	87340	N	-			
HCAB	HEPATITIS C AB TEST	86803	N	-			
HCABCF	HEP C AB TEST, CONFIRM	86804	N	-			
HCGQLS	HCG SERUM QL	84703	N	-			
HCGQT	HCG QT	84702	N	-			
HPCRQL	HEPATITIS C, RNA, AMPPR QL	87521	N	-			
HCRNAQNT	HEPATITIS C, RNA, QT W/OUT GENO. REFLEX	87522	N	-			
HCRNAQT	HEPATITIS C, RNA, QT W/ GENO. REFLEX	87522	N	-			
<i>REFLEX COMPONENT:</i>	GENOTYPE, DNA, HEPATITIS C	87902	N	-			
HCT	HEMATOCRIT	85014	N	-			
HCTSP	HEMATOCRIT, SPUN	85013	N	-			
HDL	HDL CHOLESTEROL	83718	N	-			
HEMOP	HEMOGLOBINOPATHY PNL	-	Y	-			
<i>COMPONENT:</i>	HEMOGLOBIN, CHROMA	83021	-	1			
<i>COMPONENT:</i>	SICKLE CELL TEST RED	85660	-	1			
HEPAIGG	HEPATITIS A AB TOT	86708	N	-			
HEPARIN	HEPARIN ASSAY	85520	N	-			
HEPNEUT	HEPARIN NEUTRAL	85525	N	-			
HERHEMO	HEREDITARY HEMOCHROMATOSIS	-	Y	-			
<i>COMPONENT:</i>	MOLEC ISO EXT HP EA	83891	-	1			
<i>COMPONENT:</i>	MOLEC DIAG EA ENZYME	83892	-	3			
<i>COMPONENT:</i>	MOLEC SEP GEL ELEC EA	83894	-	3			
<i>COMPONENT:</i>	MOLEC NUC AMPLI EA	83898	-	2			
<i>COMPONENT:</i>	MOLEC INTERP-RPT	83912	-	1			
HERPPCR	HERPES (HSV) TYPES 1/2 DNA PCR	-	Y	-			
<i>COMPONENT:</i>	HSV, DNA, AMP PROBE	87529	-	2			
HERPSTY1	HERPES SIMP TYPE 1	86695	N	-			
HERPSTY2	HERPES SIMP TYPE 2	86696	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
HETPH	HETEROPHILE AB SCREE	86308	N	-			
HEXPPOS1	HEXAGONAL PHOSPHOLIPID	85598	N	-			
HFP	HEPATIC FUNCTION PAN	80076	N	-			
HGB	HEMOGLOBIN	85018	N	-			
HGBA1C	HEMOGLOBIN A1C	83036	N	-			
HGBCHROM	HEMOGLOBIN, CHROMA	83021	N	-			
HGBHCTSN	HEMOGLOBIN-HCT, SPUN	-	Y	-			
<i>COMPONENT:</i>	HEMATOCRIT, SPUN	85013	-	1			
<i>COMPONENT:</i>	HEMOGLOBIN	85018	-	1			
HISTAG	HISTOPLASMA AG, URINE (CAPSULATUM)	87385	N	-			
HIVRAPID	HIV-1 ELISA	G0433	N	-			
HIV1RQT	HIV-1 RNA QNT	87536	N	-			
HIV1WB	HIV CONFIRM BY WESTERN BLOT	86689	N	-			
HIV12RS	HIV-1-HIV-2, RAPID SCREEN	86703	N	-			
HLAABCSG	HLA A, B, OR C, SING	86812	N	-			
HLB27DA	HLA B27 DIS ASSOCIATION	-	Y	-			
<i>COMPONENT:</i>	MOLEC ISO EXT EA NT	83890	-	1			
<i>COMPONENT:</i>	MOLEC DOTSLLOTBLOT EA	83893	-	6			
<i>COMPONENT:</i>	MOLEC NUC ACID PROBE	83896	-	6			
<i>COMPONENT:</i>	MOLEC NUC AMPLI EA	83898	-	1			
<i>COMPONENT:</i>	MOLEC INTERP-RPT	83912	-	1			
HMSC	HEAVY METALS, SCREEN	83015	N	-			
HOMOCY	HOMOCYSTEINE	83090	N	-			
HOMOVANA	HOMOVANILLIC ACID	83150	N	-			
HPA	HPYLORI AG, STOOL, EIA	87338	N	-			
HPYL	H PYLORI	86677	N	-			
HPYLIGG	H PYLORI IGG	86677	N	-			
HPYLIGM	H PYLORI IGM	86677	N	-			
HPYLPID	H PYLORI PRESUMP ID	87081	N	-			
HSNS	HERPES SIMP NON-SPEC	86694	N	-			
HSVULT	VIRUS ISOLATION GEN	87255	N	-			
HSVDA	HERPES SIMPLEX DNA AMP PROBE	87529	N	-			
HSVIGM	HERPES SIMP VIRUS 1&2, IGM EIA	86694					
HTLVI	HTLV-HIV CONFIRM	86689	N	-			
HTLV1AB	HTLV-I AB	86687	N	-			
HTLV2AB	HTLV-2 AB	86688	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
HYDCODON	OPIATES EACH PROCEDURE	83925	N	-			
IDOUCH	IMMUNODIFF, OUCHTER	86331	N	-			
IEPO	IMMUNOELECTRO, OTH	86325	N	-			
IEPS	IMMUNOELECTRO, SERUM	86320	N	-			
IGA	IGA	82784	N	-			
IGE	GAMMAGLOBULIN IGE	82785	N	-			
IGF	SOMATOMEDIN	84305	N	-			
IGFBP-3	INSULIN-LIKE GROWTH FACTOR BINDING PR3	83520	N	-			
IGG	GAMMAGLOBULIN IGG	82784	N	-			
IGGCSF	IGG CSF	82784	N	-			
IGGPN	IMMUNOGLOBULINES A/G/M, QN	82784	N	-			
IGGSUB	IGG SUB 1,2,3 OR 4	82787	N	-			
IGM	IGM	82784	N	-			
IMCMLPX	IMMUNE COMPLEX ASSAY	86332	N	-			
IMELECUR	IMMUNOELECTROPHOR UR	86334	N	-			
IMIP	IMIPRAMINE	80174	N	-			
IMMEL	IMMUNOFIXATION ELECT	86334	N	-			
IMMUELS	IMMUNOELECT S	-	Y	-			
<i>COMPONENT:</i>	GAMMAGLOB M-G-A-D	82784	-	3			
<i>COMPONENT:</i>	IMMUNOFIXATION ELECT	86334	-	1			
IMMUNOS	IMMUNOASSAY, NOS	83520	N	-			
IMMUNRIA	IMMUNOASSAY, RIA	83519	N	-			
INDIKPR	INDIA INK PREP	87210	N	-	LabCorp via MWMC		
INFAGIQT	INF AGENT IMMUNO QUA	86317	N	-			
INFLAAB	INFLUENZA A AB BLD	86710	N	-			
INFLUAB	INFLUENZA VIRUS AB	86710	N	-			
INFLUWOP	INFLUENZA ASSAY W-OP	87804	N	-			
INHIBINB	INHIBIN B (IMMUNOASSAY, NOS)	83520	N	-			
ININK	CULTURE MYCOBACTERIA	87116	N	-			
INSULIN	INSULIN, TOTAL	83525	N	-			
INTER	RECEPTOR, NON-ENDO	84238	N	-			
INTFAB	INTRINSIC FACTOR AB	86340	N	-			
IRON	IRON	83540	N	-			
K	POTASSIUM, SERUM	84132	N	-			
KAPLAML	KAPPA LAMBDA LIGHT	-	Y	-			
<i>COMPONENT:</i>	PHELOMETRY, NOS	83883	-	2			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
KGS	KETOGENIC STEROIDS	83582	N	-			
KOHPREP	KOH PREP	87210	N	-			
LA	LACTIC ACID	83605	N	-			
LAMOTRIG	LAMOTRIGINE	80299	N	-			
LAPWCOU	LEUKO ALK PHOS W-COU	85540	N	-			
LCMV	LYMPH CHORIOMENINGITIS	86727	N	-			
LDH	LDH	83615	N	-			
LDHISO	LDH ISOENZYMES	83625	N	-			
LDL	LDL CHOLESTEROL	83721	N	-			
LEAD	LEAD	83655	N	-			
LEGPAG	LEGION PNEUMOPH AG, I	87278	N	-			
LEGGIGG	LEGIONELLA AB	86713	N	-			
LEGGIGM	LEGION PNEUM IGM	86713	N	-			
LEVETRA	LEVETIRA-KEPPRA	80299	N	-			
LH	LEUTEINIZING HORMONE (LH)	83002	N	-			
LI	LIDOCAINE	80176	N	-			
LIPASE	LIPASE	83690	N	-			
LIPID	LIPID	80061	N	-			
LITHIUM	LITHIUM	80178	N	-			
LKMAB	LIV-KID MICRO ANTI	86376	N	-			
LPHR	LIPOPROTEIN SUBFRACT PROF VAP	-	Y	-			
<i>COMPONENT:</i>	LIPOPROTEIN, HIGH RES	83701	-	1			
<i>COMPONENT:</i>	TRIGLYCERIDES	84478	-	1			
LPPLA2	LIPOPROTEIN, PL-A2	84999	N	-			
LUPPN	COAG - LUPUS INHIBITOR PANEL	-	Y	-			
<i>COMPONENT:</i>	RUSSELL VIPER VEN DI	85613	-	1			
<i>COMPONENT:</i>	THROMBOPLAST TIMEPTT	85730	-	1			
LYMEAB	LYME AB TOTAL/IGM RESPONSES	86618	N	-			
LYMEDP	LYME DIS DNA AMP PRO	87476	N	-			
LYMEIGM	LYME AB, IGM	86618	N	-			
LYMEWBPN	LYME DISEASE, WESTERN BLOT	86617	N	-	LabCorp via MWMC		
LYMS3PN	CD4 & CD8 LYMPH & T-CELL SUBSETS	-	Y	-			
<i>COMPONENT:</i>	T CELLS, TOTAL	86359	-	1			
<i>COMPONENT:</i>	T CELL, ABSO CD4-CD8	86360	-	1			
LYRICA	CHROMATOG SING - PREGABALIN	82491	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
LYTES	ELECTROLYTE	80051	N	-			
MACXPAR	PARASITE ID, MACRO-MICRO	87169	N	-			
MAGRBC	MAGNESIUM RBC	83735	N	-			
MALBLDSM	MALARIA BLOOD SMEAR	87207	N	-			
MAVINDA	M.AVIUM-INTRA, DNA, AMP PROB	87561	N	-			
MDIFF	WBC DIFF, MAN	85007	N	-			
MERCURY	MERCURY	83825	N	-			
METANEPH	METANEPHRINES, PLASMA (FREE)	83835	N	-			
METHAD	METHADONE	83840	N	-			
METHAN	VOLATILES	84600	N	-			
METHMAL	METHYLMALONIC ACID	83921	N	-			
METHTRX	METHOTREXATE	80299	N	-			
MG	MAGNESIUM	83735	N	-			
MICAB	MICROSOMAL AB	86376	N	-			
MICALBCR	MICROAL-CREAT W-RATIO, UR	-	Y	-			
COMPONENT:	MICROALB QNT	82043	-	1			
COMPONENT:	CREATININE UR-OTHER	82570	-	1			
MITCHAB	MITOCHONDRIAL AB	86255	N	-			
MMRPN	MMR TITER	-	Y	-			
COMPONENT:	MUMPS AB	86735	-	1			
COMPONENT:	RUBELLA AB IGG	86762	-	1			
COMPONENT:	RUBEOLA AB	86765	-	1			
MN	MANGANESE	83785	N	-			
MPNEDA	MYCOPLASMA PNEUMONIAE, DNA PCR	87581	N	-			
MTBDNAMP	M. TUBERCULO, DNA, AMPP	87556	N	-			
MTHFR	MTHFR MUTATION DETECTION - C677T & A1298C	81291	N	-	LabCorp via MWMC	WHOLE BLOOD OR BUCAL SWAB	refrigerate
MUMPAB	MUMPS AB	86735	N	-			
MUMPIGG	MUMPS AB IGG EIA	86735	N	-			
MUMPIGM	MUMPS IGM AB	86735	N	-			
MXSTPTT	MIXING STUDY PTT	-	Y	-			
COMPONENT:	THROMBOPLAST TIMEPTT	85730	-	1			
COMPONENT:	PTT SUBSTITUT	85732	-	1			
MYCOAB	MYCOPLASMA AB	86738	N	-			
MYCOIGG	MYCOPLASMA IGG AB	86738	N	-			
MYCOIGM	MYCOPLASMA IGM AB	86738	N	-			
MYOQT	MYOGLOBIN	83874	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
NA	SODIUM, SER PL OR WB	84295	N	-			
NAOTH	SODIUM OTHER SOURCE	84302	N	-			
NEPHLMTR	NEPHELOMETRY, NOS	83883	N	-			
NEUROLV	NEURONTIN LEVEL	80299	N	-			
NICKEL	NICKEL	83885	N	-			
NKCELTCT	NK CELLS, TOTAL CNT	86357	N	-			
NLKVAG	NORWALK LK VIR AG	87449	N	-			
OASQT	ORGAN ACID, SIN, QT	83921	N	-			
OLIGO	OLIGOCLONAL BANDS	83916	N	-			
OSMOL	OSMOLALITY, BLOOD	83930	N	-			
OVAPN	OVA & PARASITES	-	Y	-			
<i>COMPONENT:</i>	OVA-PARASITES SMEARS	87177	-	1			
<i>COMPONENT:</i>	SMEAR COMPLEX STN	87209	-	1			
OXCRBZP	OXCARBAZEPINE	82491	N	-			
PAGGLUTS	PART AGGLUT SCREEN	86403	N	-			
PAPTEST	GYN PAP TEST, LIQUID-BASED PREP	88142	N	-			
PANCA	P ANCA	86021	N	-			
PANCAIGG	PANCA IGG	83520	N	-			
PARVPN	PARVOVIRUS-IGG, IGM	-	Y	-			
<i>COMPONENT:</i>	PARVOVIRUS AB	86747	-	2			
PARV19PC	PARVOB-19 PCR	87798	N	-			
PCALTPCT	PROCALCITONIN	84145	N	-	LabCorp via MWMC	red or green	serum or plasma, frozen
PCARAGI	PNEUMOCYSTIS CAR DFA	87281	N	-			
PCPQL	PHENCYCLIDINE SERUM	G0431	N	-			
PFA	PLATELET AGGREGATION (PFA)	85576	N	-			
PHENOB	PHENOBARBITAL	80184	N	-			
PHENY	PHENYTOIN TOTAL	80185	N	-			
PHENYFR	PHENYTOIN FREE	80186	N	-			
PHENYLAL	PHENYLALANINE, BLOOD	84030	N	-			
PHOACT	ACID PHOSPHATASE	84060	N	-			
PHOSGLYC	PHOSPHATIDYLGLYCEROL	84081	N	-			
PHOSRAUR	PHOSPHORUS RANDOM UR	84105	N	-			
PINWORM	PINWORM EXAM	87172	N	-			
PLASMDAB	PLASMODIUM AB	86750	N	-			
PLT	CBC-PLT AUTORP CDM	85049	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
PLTNOTFY	PLATELET BLOOD BANK 1 UNIT	P9019	N	-			
<i>SUPP.CHARGE:</i>	S-H PLAT PHER	P9034	N	-			
<i>SUPP.CHARGE:</i>	S-H PLAT EACH UNIT	P9019	-	-			
<i>SUPP.CHARGE:</i>	S-H PLAT PHER	P9034	-	-			
<i>SUPP.CHARGE:</i>	S-H PLAT PHERE L-R	P9035	-	-			
<i>SUPP.CHARGE:</i>	S-H PLAT PHERE IRRAD	P9036	-	-			
PLTTNEUT	PLTT NEUTRALIZATION	85597	N	-			
POCPTINR	PROTHROMBIN TIME-POC INR	85610	N	-			
PO4	PHOSPHORUS, BLOOD	84100	N	-			
PPLAABPN	PARANEOPLASTIC ANTIBODIES	-	Y	-			
<i>COMPONENT:</i>	IMMUNOASSAY, RIA	83519	N	5			
<i>COMPONENT:</i>	IMMUNOASSAY, NOS	83520	N	1			
<i>COMPONENT:</i>	FLUOROSCENT AB, TITER	86256	N	9			
PRCAC	PROTEIN C ACTIVATED	85307	N	-			
PRCRAU	PROT-CREAT RATIO UR	-	Y	-			
<i>COMPONENT:</i>	CREATININE UR-OTHER	82570	-	1			
<i>COMPONENT:</i>	PROT TOTL RANDOM URI	84156	-	1			
PREALB	PREALBUMIN	84134	N	-			
PRMIHACA	INFLIXIMAB (REMICADE) HACA - PROMETHEUS	-	Y	-	PHL via MWMC	7.5mL SST	refrigerate
<i>COMPONENT:</i>	IMMUNOASSAY, NOS	83520	-	2			
PROELC	PROTEINELECTROPHORESIS FR-QNT	84165	N	-			
PROGES	PROGESTERONE	84144	N	-			
PROINS	PROINSULIN	84206	N	-			
PROLAC	PROLACTIN	84146	N	-			
PROTCAG	PROTEIN C, ACTIVITY	85303	N	-			
PROTOTOT	T PROTEIN TOT OTHR SRC	84157	N	-			
PROTS	PROTEIN S, FREE	85306	N	-			
PROTSFU	PROTEIN S, FUNCTIONAL	85306	N	-			
PRTEL	PROTEIN TOT ELECTR-INTERP, SER	-	Y	-			
<i>COMPONENT:</i>	PROTEIN TOT SER PLWB	84155	-	1			
<i>COMPONENT:</i>	PROTEINELECTR FR-QNT	84165	-	1			
PSAFRE	PSA FREE AND TOTAL	-	Y	-			
<i>COMPONENT:</i>	PSA, TOTAL	84153	-	1			
<i>COMPONENT:</i>	PSA, FREE	84154	-	1			
PSAT	PSA, TOTAL	84153	N	-			
PT	PROTHROMBIN TIME	85610	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
PTH	PARATHORMONE PTH	83970	N	-			
PTHI	PTH INTACT	83970	N	-			
PTHRP	PTH REL PROTEIN	83519	N	-			
PTPTT	COAGP	-	Y	-			
<i>COMPONENT:</i>	PROTHROMBIN TIME	85610	-	1			
<i>COMPONENT:</i>	THROMBOPLAST TIMEPTT	85730	-	1			
PTSUBST	PROTHROMBIN, SUBSTITUTE	85611	N	-			
PTTSUB	PTT SUBSTITUT	85732	N	-			
P2Y12	PLATELET AGGREGATION (PLAVIX)	85576	N	-			
QFEV	COXIELLA BRUN QFEVER	86638	N	-			
RBC	RBC COUNT	85041	N	-			
RENALDRT	RENIN-ALDOST RATIO	-	Y	-			
<i>COMPONENT:</i>	ALDOSTERONE	82088	-	1			
<i>COMPONENT:</i>	RENIN	84244	-	1			
RENIN	RENIN	84244	N	-			
REPTILAS	REPTILASE TEST	85635	N	-			
RETICAB	RETICULIN AB	88347	N	-			
RETICM	RETIC COUNT, MANUAL	85044	N	-			
RETLDPLT	RETICULATED PLT ASSAY	85055	N	-			
RFP	RENAL FUNCTION PANEL	80069	N	-			
RFQL	RHEUMATOID FACT, QL	86430	N	-			
RH	BLOOD TYPING RH D	86901	N	-			
RHOGAMWK	RHO GAM WORKUP	-	Y	-			
<i>COMPONENT:</i>	BLOOD TYPING ABO	86900	-	1			
<i>COMPONENT:</i>	BLOOD TYPING RH D	86901	-	1			
<i>COMPONENT:</i>	ANTIBODY SCREEN, RBC	86850	-	1			
<i>COMPONENT:</i>	FETAL HGB, ROSETTE	85461	-	1			
RMSFAB2	RICKETTSIA AB - ROCKY MTN FEVER	86757	N	-			
ROTAAG	ROTAVIRUS AG, EIA	87425	N	-			
RPR	RPR QL	86592	N	-			
RSV	RESP SYNCYTIAL AG, EI	87420	N	-			
RUBEAB	RUBEOLA AB	86765	N	-			
RUBEOIG	RUBEOLA IGG	86765	N	-			
RUBEOIM	RUBEOLA IGM AB	86765	N	-			
RUBIGG	RUBELLA AB IGG	86762	N	-			
RVPPCR	RESP VIRUS PCR	-	Y	-			
<i>COMPONENT:</i>	DETECT AGT NOS, DNAAM	87798	-	10			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
RVVD	RUSSELL VIPER VEN DI	85613	N	-			
SACCR	SACCHAROMYCES CEREVISIAE PANEL	-	Y	-			
<i>COMPONENT:</i>	FUNGUS AB, NOS	86671	-	2			
SAL	SALICYLATE	80196	N	-			
SCL70	ANTISCLERODERMA-70 ABS	86235	N	-			
SEROREL	SEROTONIN RELEASE ASSAY	86022	N	-			
SEROTON	SEROTONIN	84260	N	-			
SFCRY	EXAM, SYNOVIAL FLU CR	89060	N	-			
SG	SPECIFIC GRAVITY	84315	N	-			
SHBG	SEX HORMONE BINDING GLOBULIN	84270	N	-			
SICKLE	SICKLE CELL TEST RED	85660	N	-			
SJGRSYN	SJOGRENS ABS (ANTI-SS-A/SS-B)	-	Y	-			
<i>COMPONENT:</i>	ANTI-ENA	86235	-	2			
SKINSCR	TISSUE EXAM FOR FUNGI (KOH)	87220	N	-			
AGNOSWOP	STREP PNEUM AG (NOS ASSAY WOP)	87899	N	-			
SODRSCR5	CHILD ENDANGERMENT DRUG SCREEN		Y	-	PHL via MWMC		HAS OPTION TO ADD ON OPIOID PANEL IF REQ BY DR
<i>COMPONENT:</i>	UDS CHEM ANALYZER	G0431	-	5			SEE HANDLING INSTRUCTIONS
SPNEURAG	STREP PNEUMO RAPID AG	87449	N	-			
STRADPRO	STREP A, DNA, DIR PROB	87650	N	-			
STRAOPTI	STREP A ASSAY W-OPTI	87880	N	-			
STREPAAG	STREP A AG, EIA	87430	N	-			
STRMAB	SMOOTH MUSCLE AB SCRN	86255	N	-			
SUGSQT	SUGARS, SING, QT	84378	N	-			
TACRO	TACROLIMUS	80197	N	-			
TBTSTMEA	QUANTIFERON GOLD TB TEST	86480	N	-			
TCELLTCT	T CELLS, TOTAL	86359	N	-			
TESTOFR	TESTOSTERONE, FREE	84402	N	-			
TESTOTO	TESTOSTERONE, TOTAL	84403	N	-			
TESTSHBG	BIOAVAILABLE TESTOSTERONE & SEX HORM BNG GLB	-	Y	-			
<i>COMPONENT:</i>	TESTOSTERONE, TOTAL	84403	-	1			
<i>COMPONENT:</i>	SEX HORMONE BINDING GLOBULIN	84270	-	1			
TGAB	THYROGLOBULIN ANTIBODY	-	Y	-			
<i>COMPONENT:</i>	THYROGLOBULIN	84432	-	1			
<i>COMPONENT:</i>	THYROGLOBULIN AB	86800	-	1			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
TGLTSIGA	TRANSGLUTAM TISSUE IGA	83516	N	-			
TGLTSIGG	TRANSGLUTAM TISSUE IGG	83516	N	-			
THEO	THEOPHYLLINE	80198	N	-			
THROMEV	COAG STUDY THROM EVALUAYION	-	Y	-			
COMPONENT:	MOLEC ISO EXT HP EA	83891	-	1			
COMPONENT:	MOLEC NUC ACID PROBE	83896	-	2			
COMPONENT:	MOLEC NUC AMPLI EA	83898	-	1			
COMPONENT:	MOLEC MUTATION SCAN	83903	-	1			
COMPONENT:	MOLEC INTERP-RPT	83912	-	1			
COMPONENT:	BETA-2 GLYCOPROT AB	86146	-	2			
COMPONENT:	CARDIOLIPIN AB, IGG	86147	-	2			
COMPONENT:	ANTITHROMBIN IIIACTIVITY	85300	-	1			
COMPONENT:	PROTEIN C, ACTIVITY	85303	-	1			
COMPONENT:	PROTEIN S, FREE	85306	-	1			
COMPONENT:	PROTEIN C, ACTIVATED	85307	-	1			
COMPONENT:	RUSSELL VIPER VEN DI	85613	-	1			
COMPONENT:	THROMBOPLAST TIMEPTT	85730	-	1			
THYROGAB	THYROGLOBULIN AB	86800	N	-			
THYROGQT	THYROGLOBULIN	84432	N	-			
TIBC	IRON BINDING CAPACITY	83550	N	-			
TISTGG	TISSUE TRANSG IGG (IMMUNOASSAY, NOS)	83520	N	-			
TOBRAPK	TOBRAMYCIN PEAK	80200	N	-			
TOBRAR	TOBRAMYCIN RANDOM	80200	N	-			
TOBRATR	TOBRAMYCIN TROUGH	80200	N	-			
TOPR	TOPIRAMATE	80201	N	-			
TOTT3	T3, TOTAL	84480	N	-			
TOXIGG	TOXOPLASMA AB, IGG	86777	N	-			
TOXIGM	TOXOPLASMA AB, IGM	86778	N	-			
TPAB	THYROID PEROXIDASE AB	86372	N	-	LabCorp via MWMC		
TPALIDAB	TREPONEMA PALLIDUM AB	86780	N	-			
TPROTQT	PROTEIN TOTAL SERUM	84155	N	-			
TRANSF	TRANSFERRIN	84466	N	-			
TRANSFRX	TRANSFUSION REACTION	-	Y	-			
COMPONENT:	BLOOD TYPING ABO	86900	-	2			
COMPONENT:	BLOOD TYPING RH D	86901	-	2			
COMPONENT:	DIRECT COOMBS TEST	86880	-	2			
TRAZ	TRAZADONE (DRUG QUANT)	80299	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
TRIG	TRIGLYCERIDE	84478	N	-			
TROPT	TROPONIN, QT	84484	N	-			
TRSNG	TRYPSINOGEN (TRYPSIN) - IMMUNOASSAY, RIA	83519	N	-			
TRYP	TRYPTASE	83520	N	-			
TSH	TSH (THYROID STIMULATING HORMONE)	84443	N	-			
TSHIGAB	THYROID STIM IMMUNOGLOBULIN	84445	N	-			
TT	THROMBIN TIME	85670	N	-			
TUMAGO	TUMOR AG, OTHER (NEURON SPECIFIC ENOLASE)	86316	N	-			
TYPESCRN	TYPE AND SCREEN	-	Y	-			
<i>COMPONENT:</i>	BLOOD TYPING ABO	86900	-	1			
<i>COMPONENT:</i>	BLOOD TYPING RH D	86901	-	1			
<i>COMPONENT:</i>	ANTIBODY SCREEN, RBC	86850	-	1			
T3REV	T3, REVERSE	84482	N	-			
T4TOT	THYROXINE T4, TOT	84436	N	-			
T3U	T3, UPTAKE	84479	N	-			
UA	UA, AUTO, W-O SCOPE	81003	N	-			
UALB	ALBUMIN URINE	82042	N	-	Chemistry	Red or Green	No SST or PST tubes
UAM	UA NONAUTO W-O SCOPE	81002	N	-			
UAMIC	MICROSCOPIC EX OF UR	81015	N	-			
UARSENIC	ARSENIC URINE	82175	N	-			
UAWMIC	URINALYSIS W-MICRO	81001	N	-			
UCA	CALCIUM URINE RANDOM	82340	N	-			
<i>COMPONENT:</i>	CADMIUM	82300	N	-			
<i>COMPONENT:</i>	CREATININE UR-OTHER	82570	-	-			
UCATF	CATECHOL FR-U RANDOM	82384	N	-			
UCHLAMGC	UR CHLAMYDIA & GONORRHOEAE, TMA	-	Y	-			
<i>COMPONENT:</i>	CHLMYD TRACH DNA AMP	87491	-	1			
<i>COMPONENT:</i>	N. GONORRHO, DNA, AMPPR	87591	-	1			
UCL	CHLORIDE URINE RAN	82436	N	-			
UCREA	CREATININE URINE RANDOM	82570	N	-			
UDSSINGL	DRUG SCREEN SINGLE	G0431	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
UHMSC	HEAVY METAL URINE	-	Y	-			
<i>COMPONENT:</i>	CREATININE UR-OTHER	82570	-	1			
<i>COMPONENT:</i>	LEAD	83655	-	1			
<i>COMPONENT:</i>	MERCURY	83825	-	1			
<i>COMPONENT:</i>	ARSENIC	82175	-	1			
UIMMEPHO	IMM E-PHORS-URIN-CSF	86335	N	-			
UK	POTASSIUM, URINE	84133	N	-			
ULEGAG	LEGIONELLA AG UR	87449	N	-			
UMETA	METANEPHRINS URINE	83835	N	-			
UMYOQT	MYOGLOBIN URINE	83874	N	-			
UNA	SODIUM, URINE RANDOM	84300	N	-			
UOSMOL	OSMOLATITY, URINE	83935	N	-			
UOXAL	OXALATE URINE, 24 HR	83945	N	-			
UPORQT	PORPHYRINS, UR, QT	84120	N	-			
UPREG	HCG QL URINE	81025	N	-			
UPRO	PROTEIN TOTL RANDOM UR	84156	N	-			
UPROTPHO	PROTEIN ELEC URIN CSF	84166	N	-			
URIC	URIC ACID, BLOOD	84550	N	-			
UURIC	URIC ACID, URINE RAND	84560	N	-			
UVMA	VMA, URINE	84585	N	-			
UVOL	URINALYSIS, VOL MEASU	81050	N	-			
U24CA	CALCIUM UR 24HR	82340	N	-			
U24CAD	CADMIUM, UR 24HR OR RANDOM	82300	N	-	LabCorp via MWMC	PLASTIC CONTAINER	room temp
U24CIT	CITRIC ACID, 24 HR URINE	82507	N	-			
U24CL	CHLORIDE URINE 24H	82436	N	-			
U24CORTF	CORTISOL FREE UR 24HR	82530	N	-			
U24CRE	CREATININE UR 24HR	82570	N	-			
U24K	POTASSIUM UR 24HR	84133	N	-			
U24LYTES	ELECTROLYTE URINE 24	-	Y	-			
<i>COMPONENT:</i>	POTASSIUM, URINE	84133	-	1			
<i>COMPONENT:</i>	SODIUM, URINE RANDOM	84300	-	1			
<i>COMPONENT:</i>	CHLORIDE URINE RAN	82436	-	1			
U24MERC	MERCURY URINE, 24HR	83825	N	-			
U24MET	METANEPHRINS URINE, 24HR	83835	N	-			
U24MICAL	MICROABLUMIN, 24HR URINE	82043	N	-			
U24NA	SODIUM UR 24HR	84300	N	-			
U24NITRO	UREA NITROGEN UR RDM OR 24HR	84540	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
U24PBGQT	PORPHOBILINOGEN, QT UR 24HR	84110	N	-			
U24PHOS	PHOSPHORUS 24HR UR	84105	N	-			
U24PRO	URINE PROTEIN UR 24HR	84156	N	-			
U24URIC	URIC ACID URINE 24	84560	N	-			
VALP	VALPROIC ACID	80164	N	-			
VANPK	VANCOMYCIN PEAK	80202	N	-			
VANR	VANCOMYCIN RANDOM	80202	N	-			
VANTR	VANCOMYCIN TROUGH	80202	N	-			
VARIDNA	VARICELLA DNA PCR	87798	N	-			
VARZAB	VARICELLA-ZOSTER AB	86787	N	-			
VARZAGI	VARICELLA ZOSTER AGI	87290	N	-			
VARZG	VARICELA ZOSTR IGG	86787	N	-			
VARZM	VARICELA ZOSTER IGM	86787	N	-			
VASOPR	VASOPRESSIN	84588	N	-			
VIP	VASOINTESTINAL PEP	84586	N	-			
VIRABNO	VIRUS AB, NOS	86790	N	-			
VISC	VISCOSITY EXAM	85810	N	-			
VITA	VITAMIN A	84590	N	-			
VITB1	VITAMIN B-1 (THIAMINE), PLASMA	84425	N	-			
VITB1WB	VITAMIN B-1 (THIAMINE), WHOLE BLOOD	84425	N	-			
VITB12	VITAMIN B-12	82607	N	-			
VITB12FA	VIT B 12-FOLIC ACID	-	Y	-			
<i>COMPONENT:</i>	VITAMIN B-12	82607	-	1			
<i>COMPONENT:</i>	FOLIC ACID SERUM	82746	-	1			
VITB6	VITAMIN B-6	84207	N	-			
VITC	VITAMIN C	82180	N	-			
VITD325O	VITAMIN D-3 25-OH	82306	N	-			
VITE	VITAMIN E	84446	N	-			
VITK	VITAMIN K	84597	N	-			
VONWIL	FACTOR VIII VW AG	85246	N	-			
WBC	WBC COUNT	85048	N	-			
WBCCTBF	WBC CT BODY FLUID	89051	N	-			
WBCFECAL	WBC STOOL (FECAL LEUKOCYTES)	89055	N	-			
WETPREP	SMEARWET MOUNT SAL-I	87210	N	-			
WNILEAB	WEST NILE VIR AB	86789	N	-			
WNILEIGG	WEST NILE VIR AB IGM	86788	N	-			

TEST CODE	TEST	CPT	PANEL	QTY	PROCESSING SITE	TUBE TYPE	ADDITIONAL INFO
WNILRPCS	W. NILE RNA-PCR CSF	87798	N	-			
ZINC	ZINC	84630	N	-			
ZINCP	PROTOPORPH, RBC, QT	84202	N	-			
17KETO	17-KETOSTEROIDS, TOT	83586					
17OHPROG	HYDROXYPROGES 17-D	83498	N	-			
5HIAA	5-HIAA	83497	N	-	PHL	Sterile Container	80455 Send aliquot - note total volume on order
6TGNMMPN	THIOPURINE METABLOTIES (TPMT) - PROMETHEUS	82491	N	-	PHL		
BONEMARR	BONE MARROW NOTICE						
CYTORDER	CYTOLOGY ORDER						
PTHORDER	PATHOLOGY ORDER						